Universal composite

INSTRUCTIONS FOR USE

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With non-stick handling, the radiopaque light cured, condensable and low-shrinkage composite, LUNA 2 offers polishability and strength on an unprecedented scale for anterior and posterior restorations alike.

COMPOSITION:

24% wt (44% vol.) multifunctional methacrylic ester,
76% wt (56% vol.) inorganic filler (40nm - 1.5 micron).

INDICATIONS:

- Direct anterior/posterior restorations and veneers.
- Core build ups.
- Indirect inlays, onlays and veneers.
- Splinting.
- Composite and porcelain repair.
- Sandwich technique with glass ionomers.

CONTRAINDICATIONS:

Pulp capping. Product may cause skin irritations to some people. In such cases, discontinue use and seek medical attention. Do not use in conjunction with any eugenol-containing materials. Patients with known sensitivity to acrylics or allergic reactions to the material of this product should not use it.

INSTRUCTIONS FOR USE:

1. Select composite shade by matching with a clean and moist tooth using VITA classical A1-D4® composite shade guide.
2. Clean and isolate tooth. Rubber dam is the preferred mode of isolation.
3. Prepare a conservative cavity, employing standard techniques.
4. Use a calcium hydroxide liner and appropriate glass ionomer cement liner to protect deep cavity preparations in close proximity to the pulp.
5. Direct adhesive application: Follow the technique recommended by the adhesive you are using (e.g. Zipbond Universal, SDI Limited)

Self-Etch Technique

Bottle delivery: Dispense one or two drops (depending on restoration size) of self-etch adhesive, such as Zipbond Universal (SDI Limited) into a mixing well. Immediately replace the cap.

Unit dose delivery: Open the unit dose device and insert a brush applicator (e.g. SDI Points) into the container to saturate the tip.

NOTE: Proceed to the next steps quickly as the liquid can evaporate and result in poor bonding.

- a) Apply adhesive into the cavity by scrubbing for 10 seconds with a brush (e.g. SDI Points).
- b) Leave on cavity surface for 10 seconds.
- c) Blow with water/oil-free air until you see no movement of the bond (minimum 5 seconds).

NOTE: A vacuum aspirator can be used to prevent adhesive liquid from scattering.

- d) Light cure for 10 seconds (440-480nm wavelength, high power LED curing light such as Radii Xpert - SDI Limited).

NOTE: Enamel subjected to fluoridation should be etched for 90 to 120 seconds.

IMPORTANT: If using an etch syringe, ensure that the dispensing tip hub is firmly attached to the syringe by twisting the hub securely onto the syringe. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the etchant's viscosity and rate of extrusion.

Selective Etch Technique

- a) Apply phosphoric acid etching gel (e.g. SDI Super Etch 37%) to tooth enamel only and leave for 10 seconds.
- b) Rinse thoroughly with water and dry with water/oilfree air or with cotton pellets. Do not over dry.
- c) Apply adhesive, such as Zipbond Universal (SDI Limited) according to steps in "Self-Etch Technique".

Total Etch Technique

- a) Apply phosphoric acid etching gel 37% (e.g. SDI Super Etch) to tooth structure (enamel & dentine) and leave for 10 seconds.
 - b) Rinse thoroughly with water and dry with water/oilfree air or with cotton pellets. Do not over dry.
 - c) Apply adhesive according to "Self Etch Technique".
6. Place composite in increments of 2mm or less.

NOTE: If using the syringe system, after each use, turn the screw plunger anti-clockwise half a turn to release any built up pressure.

7. Light-cure each increment for at least 20 seconds for light shades (all shades except OA2, OA3, OA3.5) and at least 40 seconds for opaque shades (OA2, OA3, OA3.5).
8. Shape, polish and finish employing standard techniques.

OPTION: Additional polishing using polishing paste:

After step 8, polishing paste may be used to give the composite a lustrous finish, using the following steps:

- a) Apply Polishing Paste (SDI Limited) to a rubber cup or disc.
- b) Smear a thin layer of paste onto the restoration.
- c) Add a small amount of water to both the tooth and cup. Polish for 30 seconds at low speed and light pressure.
- d) Wash and dry.

STORAGE AND HANDLING:

- Store at temperatures between 10° and 25°C (50° - 77°F).
- Use at room temperature between 20° and 25°C (68° - 77°F).
- Syringe system: when not in use replace the cap tightly.

PRECAUTIONS:

- Avoid prolonged contact of etchant, adhesive and composite with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization.
- Patients with known material sensitivity or allergic reactions to this product should not use it.
- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID:

- **Eye (contact):** Wash thoroughly with water and seek medical attention.
- **Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- **Ingestion:** Drink plenty of water/milk. Seek medical attention if symptoms persist.
- **Inhalation:** No symptoms expected.

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SDI Holdings PTY Ltd do Brasil
Made in Australia by SDI Limited
SDI Germany GmbH